## TIRAS DE GLICOSE / TIRAS DE GLUCOSA / GLUCOSE STRIPS

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## ENGLISH

### USE

Screening test for early diagnosis of diabetes, suitable for serial examination.

### PRINCIPLE

The test is based on a double sequential enzyme reaction. One enzyme glucose oxidase catalyzes the formation of gluconic acid and hydrogen peroxide from the oxidation of glucose. A second enzyme peroxidase catalyzes the reaction of hydrogen peroxide with potassium iodide chromogen to oxidize the chromogen to colours ranging from blue/green to brown.

### REACTIVE INGREDIENTS (based on dry weight at time of impregnation of 100 strips)

Glicose oxidase	1,62 mg
Peroxidase	0,65 mg
Potassium iodide	9,10 mg

### PRECAUTIONS

For diagnostic use only.

Only use well washed and clean vessels for urine collection. The presence of usual urine preservatives will not affect the test results. Instructions of use should be followed correctly in order to obtain reliable test results.

The effect of medicaments or their metabolic products on the test is not known in all cases. In case of doubt it is recommended not to take the medicaments and then repeat the test.

Remove only as many test strips as are required, and reseal the container immediately after use. Do not touch the test paper. Avoid exposing the strips to sunlight and moisture.

## STORAGE

Store the container below +30°C in a dry place. The test strips are stable, when stored properly up to the date of expiry indicated.

## COLLECTION AND PREPARATION OF THE SAMPLE

Use fresh collected urine (not more than 4 hours) in a clean and well washed container .

## DIRECTIONS (FOR PATIENTS)

If test results appear doubtful please consult your doctor. This, obviously, is also required even if test results are negative.

## INSTRUCTIONS FOR USE

The test should be performed about 2 hours after breakfast.

Remove test strips from the container just before the test, and immediately reclose box tightly. Do not touch test paper.

Dip the reagent strip for approximately 2 seconds into the fresh urine. Draw it across the rim of the container to remove excess urine. After 30 to 60 seconds compare the test strip with the colour scale. The best time for comparison is after 30 seconds. Colour changes take place after more than 2 minutes are of no significance. The urine should not be more than 4 hours old when tested.

## EVALUATION

If the test field for glucose remains unchanged, the urine does not contain measurable concentrations of glucose.

Pathological glucose concentrations are indicated by a colour change blue/green to brown. The colour fields correspond to the following ranges of glucose concentrations:

0 (neg.), 100 ( $\pm$ ), 250 (+), 500 (++) , 1000 (+++), 2000 (++++) mg/dl or

0 (neg.), 6 ( $\pm$ ), 15 (+), 30 (++) , 60 (+++), 120 (++++) mmol/L

## QUALITY CONTROL

For best results, reagent strips performance should be confirmed with samples and/or positive and negative controls, every time the new container is opened.

## LIMITATIONS

In any case, in order to establish a final diagnosis and prescribe the appropriate therapy, your doctor should verify results obtained with test strips, with diagnostic methods of his own.

## SOURCES OF ERROR

Large amounts of ascorbic acid may be present in urine after a high intake of vitamin C (e.g. vitamin tablets, antibiotics or fruit juices) can lead to lower or falsely negative results. In addition an inhibitory effect is produced by gentisic acid. Falsely positive reactions can also be produced by a residue of peroxide containing cleaning agents.

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